

Verbund für Angewandte Hygiene (VAH)  
Association for Applied Hygiene

- Desinfektionsmittel-Kommission im VAH -  
Disinfectant Commission in VAH

## Zertifikat / Certificate

über die Konformität der Wirksamkeitsprüfungen für / for conformity of efficacy tests for

### Ster 1 plus

mit dem Anforderungskatalog der Deutschen Gesellschaft für Hygiene und Mikrobiologie (DGHM) - Stand 4.2.2002 - basierend auf den Standardmethoden der DGHM - Stand 1.9.2001 - bzw. den Übergangsbestimmungen vom 4.9.2002. DGHM und VAH bestätigen die Äquivalenz mit den früher herausgegebenen DGHM-Zertifikaten.  
with the Requirements issued by the German Society for Hygiene and Microbiology (DGHM) on February 4<sup>th</sup> 2002 which is based on the DGHM Standard Methods of September 1<sup>st</sup> 2001 or the transitional provisions of Sept. 4<sup>th</sup> 2002. DGHM and VAH confirm the equivalence to the former issued DGHM certificates.

ANTRAGSTELLER / APPLICANT:

**CEFLA S.C.**  
**Via Selice Provinciale 23/A**  
**IT-40026 Imola (BO)**

WIRKSTOFFE nach Art und Menge bezogen auf 100 g:

Quantity of active substances per 100 g:

**25,0 g n-Propanol**

**15,0 g Ethanol**

**1,5 g Bardac 114®**


Hiermit wird bestätigt, dass das o.g. Produkt für die prophylaktische Desinfektion in den aufgeführten Anwendungsbereichen in folgenden Konzentrations-Zeit-Relationen als wirksam eingestuft wird:  
This is to confirm that the above product was found to be effective for prophylactic disinfection in the application domains listed below at the specified concentration/contact time ratios:

Flächendesinfektion zur Prophylaxe in Krankenhaus und Praxis (Bakterizidie und Levurozidie) Surface disinfection for prophylaxis in hospital and primary healthcare (bactericidal activity and yeasticidal activity)									
Wischen Mechanical action		Organische Belastung Interfering substance		Einwirkzeit in min Contact time in min					
mit/with	ohne/ without	gering/clean conditions	hoch/dirty conditions	5	15	30	60	240	
X			X	konz.					
Zusätzlich wirksam gegen / Additionally active against Mycobacterium terrae (Tuberkulozidie / tuberculocidal activity)									
X			X	konz.					
Zusätzlich wirksam gegen / Additionally active against Mycobacterium avium (Mykobakterizidie / mycobactericidal activity)									
X			X	konz.					

Das Zertifikat ist gültig vom 29.07.2014 bis zum 29.07.2017 (3 Jahre) /  
Certificate is valid from 29.07.2014 until 29.07.2017 (3 years)

Der Antragsteller hat sich mit den Bedingungen der zum Zeitpunkt der Antragstellung gültigen Geschäftsordnung der Desinfektionsmittelkommission im VAH einverstanden erklärt und rechtsverbindlich bestätigt, dass das von ihm in Handel gebrachte Präparat in seiner Zusammensetzung identisch mit den für die Erstellung der Gutachten eingereichten Mustern ist.  
The applicant has agreed to the conditions laid down in the rules of the disinfectant commission in the VAH valid at the time of application and has legally binding confirmed that the distributed product is identical with the product used for the activity testing.

Bonn, den **29.07.2014**  
Place/Date

  
Der Vorsitzende der / The Chairman of the  
Desinfektionsmittel-Kommission im VAH / Disinfectant Commission in VAH - Ident.No N14/142